

Baffin AUTOMATIC

Advanced set up guide

1. Before transferring the user in, get the pelvic belts ready laying flat on the seat.



2. Place the user in the seating position (slight tilt may be used). Ensure that the back is in the right position (good contact of the lower back and the backrest).
3. Going round the hip, fasten the pelvic belts to the black strap.



4. Fasten the vest (harness) to the coloured (usually the same colour as the equipments upholstery) strap. (Top straps are coming from the top of the backrest, bottom corner straps are coming from under the seat and should go round the hip)



5. Fasten the abduction belts if needed.



6. Fasten the straps on the foot platforms.



Now the user is secure in the seating position.

7. If we want to get the user standing, make sure all the steps before are done. When the user is seating secure, press the top left button on the remote and make sure the user is sitting in 90-90 degrees position and the seat is parallel to the floor before going to lay down.



8. When the user is sitting in 90-90 degrees position with the seat parallel to the floor, press the bottom right button to get the user to lay down. Sometimes you might need to tighten the vest (harness) and pelvic belts straps.




9. If the knees are a bit flexed get the footplates lower.



10. Get the knee blocks in. Pull the silver pin down, push the knee blocks in, let go the silver pin and wait for a click to make sure that it went securely in. Don't push the knee blocks to hard, make sure they're not too tight just below the kneecap.



11. The user is secure to stand up. Press the top left button on the remote and get user in the standing position. Don't go over 90 degrees (we recommend going to maximum around 80 degrees). Put the tray in for a better posture.



12. If you would like to go back to seating position, press the top right button to lay down, take the knee blocks out (the equipment won't let you go to seating position with the knee blocks in for safety). Then press bottom left button to go to seating position. If you would like to use tilt in space press top right button on the remote.